soulima.®

USER MANUAL

Handheld wireless massager Soulima 21630

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00021630 M: YT-1706AL

EN

Ladies and Gentlemen, thank you for purchasing our product!
Before starting to use, please read the following instructions
for the correct use of the product. Please keep these instructions on
future and to follow its recommendations, as non-compliance may endanger life or health.

Important: KEEP FOR FUTURE REFERENCE!

- Please read it completely before use and follow all instructions in the manual.
- Do not use the product if any part is missing or damaged.
- Please keep this manual for future reference.

BASIC INFORMATION:

- To reduce the risk of personal injury, electric shock, fire, equipment damage, or malfunction, read all important information and warnings in this manual before use.
- It is strictly forbidden to use this product by injured, sick and minors. If you suffer from a muscle injury, consult your doctor before using the device.
- Do not use this device on wet or dirty skin. The gun can be used in areas covered by clothing. The massager should be pressed against the body and gently moved between different points, stopping in each place for about 60 seconds.
- Do not use the massager on parts of the body above the neck or on areas with less muscle or bone content. The device may only be used on soft tissues if the user does not experience pain or other adverse reactions.
- Do not use the device on parts of the body with bruises or if there is pain or unusual body reactions.
- Do not submerge the product in water or allow water to enter the ventilation openings of the product.
- Do not place the device in liquids. Do not allow any liquid to enter the ventilation openings.
- Keep away from fire and sources of heat.
- Do not throw the device and avoid accidental falls.
- Do not disassemble or modify the product.
- Fingers, hair or other body parts may become trapped in the shaft and back of the massage head. Do not come in close contact with these areas.
- The device must not be used by children under 8 years of age. Children must not use, clean and maintain the device without adult supervision.
- Remember that if you feel uncomfortable, stop using it immediately.

People who meet any of the following conditions should consult a doctor before using the massager:

- Pregnant women
- People with a pacemaker
- People who have recently had surgery
- People with spinal, joint or bone defects
- People who have complications caused by other conditions (such as neuropathy or damage to the retina).

Modern scientific research focuses on determining the impact of vibration massage on people with special health conditions and may limit or extend the scope of the above-mentioned contraindications. Practical experience has shown that the combination of treatment and vibration massage can also serve as a type of therapy in many circumstances. However, such an approach should be defined and ordered by a medical specialist or therapist, or conducted under special supervision in their presence.

PRECAUTIONS:

- It is forbidden to run this product without the massage head installed.
- Before use, please check whether the product can work normally according to the operation method of this manual.
- Before each use, carefully check that the appliance or tip is undamaged and can be used normally according to the methods specified in the instruction manual.
- Please use this product in strict accordance with the massage position and operating method described in this manual.
- The room temperature of the usage environment is 0-40°C.
- The room temperature of the charging environment is 0-40°C.
- The time of using the same location should not exceed 60s.
- Do not use it continuously for more than 1 hour.
- Clean the device after use to avoid the accumulation of residues such as sweat or grease.
- Always disconnect the charger when charging is complete.
- Do not use the product while charging.
- Always follow the instructions for use.

SPECIFICATION:

Type: Hand massager
Number of operating modes: 12
Number of power levels: 10
Power: 24W
Voltage: 9V
Power supply: battery
Working time: 1.5-2h

Included:

massager
6 interchangeable tips
Charger

CHARACTERISTIC:

A hand-held massager with interchangeable heads is primarily used to regenerate muscles, reduce stress and tension, and relax. It will help soothe the painful effects of strength training, help to massage frozen and tense muscles. It will also work perfectly as a relaxation device. The massager is an optimal complement to the training and relaxation cycle for people who exercise regularly. People with less physical activity can use it to effectively stimulate muscle work.

The massager is designed to massage the whole body and stimulate blood flow. It has a strengthening effect on the muscles of the back and body. Massage helps to restore the proper function of nerve connections, while improving metabolism. It improves general well-being and relieves stress. Thanks to the use of heads with protrusions, the massager also serves as an effective tool in the fight against cellulite. Improves firmness, vascularity and overall appearance of the skin.

DEVICE DESCRIPTION

1. Replaceable massage heads
2. Display showing battery level and operating mode
3. Massage mode selection
4. Increasing the power of the massage
5. Reducing the massage power
6. On/Off
7. Charging

ADVANTAGES OF THE DEVICE

Massage of all parts of the body,

Capacious battery for wireless operation

Adjustment of the operating mode and massage intensity:

- 12 massage programs,
- 10 intensity levels,

6 multifunctional massage heads,

Ergonomic shape, facilitating massage in hard-to-reach places

non-slip handle,

Muscle regeneration and relaxation

APPLICATION OF MASSAGE HEADS

1. Smooth: it is used for a deep massage (shiatsu), which can be used on the neck, shoulders, back and in the area of the greatest muscle tension.

2. With an inset: an in-depth shiatsu massage allows you to relax and regenerate tense muscles even better.

3. Multi-point: ideal for deep massage of large areas of the body, such as the back, abdomen or buttocks. The head used on the thighs and buttocks also has an anti-cellulite effect.

4. Five-point: massage head that relaxes the loins and back perfectly.

5. Flat: works well when relaxing tight muscles in the thighs and calves, helps reduce swelling and improves circulation. Also suitable for shoulder massage.

6. For acupressure: a special head designed for foot and hand massage. It precisely presses small points on the body, bringing relaxation and deep relaxation. Applied to the area of thighs and buttocks, it has an anti-cellulite effect.

INSTRUCTIONS FOR USE

Landing

- Fully charge the battery before first use.
- To charge, plug the charger into a wall socket, then plug the end into the port
- charging the battery.
- When the battery is connected, the LED on the charger will turn red.
- A full charge is indicated when the LED on the charger turns green.
- The battery can be charged at any battery level.
- It is not recommended to fully discharge the battery.

Activation

Press the power button to turn on the massager.

Exclusion

To turn off the massager, press the power button.

Stages of use

Prerequisites: Before use, prepare the skin or thin clothing on the part of the muscle to be relaxed. The massaged surface should be dry and clean.

HOW TO USE:

Step 1

Install the massage head

According to the different parts of use and effects, please choose the appropriate massage head. When installing the massage head, please insert the massage head directly into the outlet of the massage gun, and insert it firmly all the way so that it aligns with the line of the massage head.

Step 2

Use a massage gun

Turn on the gun, switch to the appropriate gear and adjust the down pressure. Gently press the massage head against the muscle group to be relaxed, move it slowly back and forth, then stay there for 30 seconds after finding the sensitive and sore spot. The entire massage process should not exceed 10 minutes. The console gun itself has its own downforce, so there is no need to apply too much downforce. Avoid bones, joints, arteries and other parts.

MAINTENANCE AND CLEANING

- Pull out the massage head directly, wipe the gun head and body, clean the body with a slightly damp soft towel, and then use a cloth to dry.
- Store the product in a dry, well-ventilated and non-corrosive place.
- Do not use corrosive solutions to wipe the case, otherwise it may cause product failure or product cracking, discoloration.
- This product is suitable for steam, train, ship and airplane transportation. During transportation, the product should be fixed and placed gently to prevent violent vibration, collision, tipping over and dropping.
- When storing or moving, please clean it and put it back into the packing box or storage box.

BATTERY INFORMATION

- Keep batteries/rechargeable batteries away from children.
- Do not store batteries/rechargeable batteries in easily accessible places as there is a risk that they may be swallowed by children or pets. If swallowed, contact a doctor as soon as possible.
- Leaking or damaged batteries/rechargeable batteries may cause burns when in contact with the skin, so always wear suitable protective gloves when touching them.
- Batteries/rechargeable batteries must not be short-circuited, disassembled or thrown into fire. There is a risk of explosion!
- Do not recharge non-rechargeable batteries as there is a risk of explosion! Only recharge batteries intended for this purpose using the appropriate charger.
- Charge the batteries before inserting them into the product
- Batteries should only be charged under adult supervision.
- Do not mix alkaline, standard (carbon-zinc) and rechargeable (nickel-cadmium) batteries.
- Never mix new and old batteries.
- It is forbidden to interfere with the power cables - it may cause a short circuit.
- Exhausted batteries must be removed from the product.
- Do not dispose of used batteries in unsorted waste!
- Batteries must be disposed of in accordance with hazardous electronic waste disposal.
- If the batteries are not properly disposed of, hazardous substances may pose a risk to human health or the environment. Recycling of materials contributes to the conservation of natural resources.
- Batteries/rechargeable batteries must always be inserted with the correct polarity (plus/+ and minus/-).

Information for the consumer on the handling of waste electrical and electronic equipment:

Waste electrical and electronic equipment may not be placed with other household waste. Such equipment should be collected selectively, which is indicated by the attached marking - a crossed-out wheeled waste container. Failure to comply with the above rule, in the case of improper disposal of used equipment, may pose a threat to the natural environment and human health, resulting from the presence of dangerous components in the equipment, such as: electrical wiring, plastic, switches, etc. To avoid the above threat, these components should be collected and properly processed by specialized companies. Self-disassembly of worn equipment is unacceptable. Households play an important role in contributing to reuse and recovery, including recycling of waste equipment. This is done in particular through active participation in the waste equipment collection system.

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