Ruhhy®

USER MANUAL

Bathroom scale Ruhhy 24638

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


Watermarkly

EN

Ladies and Gentlemen, thank you for purchasing our product!
Before starting to use, please read the following instructions
for the correct use of the product. Please keep these instructions on
future and to follow its recommendations, as non-compliance may endanger life or health.

Important: KEEP FOR FUTURE REFERENCE!

- Please read it completely before use and follow all instructions in the manual.
- Do not use the product if any part is missing or damaged.
- Please keep this manual for future reference.

SPECIFICATION:

Size: 310*310*28mm,

Display screen: 95*57mm,

Units: kg/st/lb,

Load: 5kg - 225kg

Material: ABS + tempered glass + electrode sheet

CHARACTERISTIC:

The LCD analytical bathroom scale is a modern device that measures body weight, ranging from 5 kg to 225 kg. In addition to the basic weight measurement, the scale offers body composition analysis, including measurement of body fat percentage, muscle mass and water content in the body. Equipped with a clear LCD display, it also allows you to calculate the BMI (Body Mass Index). This advanced tool allows you to monitor various health parameters, supporting users in maintaining a healthy lifestyle.

BUTTON DESCRIPTION:

1. BMI reading
2. Read % body fat
3. Body water % reading
4. Read the % of muscle mass
5. Basal metabolic rate
6. Athlete mode icon
7. Man icon
8. Woman icon
9. Height displayed in feet/inches (only for lbs readout)
10. Weight reading in kg
11. Weight reading in pounds
12. Weight reading in the street
13. "Up" button
14. SET button
15. "Down" button

- "SET" button:

- Turns on the scale
- Selection confirmation button


"Up" or "Down" button

- Changes the value of height, age and activity level
- Selects gender: male/female/athlete
- Selects user P1-P10 in memory mode

PREPARATION BEFORE USE

1. Insert the battery

This scale operates on 3XAAA zinc-carbon batteries (included). Open the battery cover on the back of the scale and place the batteries into the battery compartment according to the "+" and "-" signs at the bottom of the compartment. Close the battery compartment cover. The scale will beep. If you do not intend to use this device for an extended period of time, it is advisable to remove the battery.

2. Set the units of measurement

To change the unit of measurement, press the UNIT button on the back of the scale (kg/lb/st).

3. Some models have an adhesive film on the LCD lens to prevent scratching. Remove it before use.

4. Place the scale in the correct position

Always use the scale on a flat and hard surface. Do not use on carpets. A carpet or uneven floor may affect the measurement result.

5. Start the scale


Press the scale platform to activate it. The display shows "0.0 kg" and then goes off. The scale is ready for use.

USER PROGRAMMING:

1. Press the "SET" button to turn on the scale.

2. Select memory (the scale can remember the data of 10 different people)

3. Press the "Up" or "Down" button to select the memory number (P1-10). Press the "SET" button to select the displayed memory number.

4. Select Gender or Athlete mode. Press the "Up" or "Down" button to switch between male, female or athlete. Press the SET button to confirm.

Definition of an athlete:

Trains a minimum of three times a week to improve specific skills required for a specific activity.

5. Select your activity level:

Press the "Up" or "Down" button to select your activity level (1-5). Press the SET button to select the displayed activity level.

Select your activity level according to the guidelines below:

Level 1. Sedentary/very inactive: little or no exercise

Level 2 Limited activity: exercise 1-3 days a week

Level 3 Moderate activity: exercise 4-5 days a week

Level 4 Very active: exercise 6-7 days a week

Level 5 Extremely active: physically demanding exercise or athletic training.

6. Enter your height (scale from 100 cm to 220 cm)

Press the Up or Down button to increase/decrease the height digits. Press the SET button to select the displayed height.

7. Enter your age (scale from 18 to 80)

Press the Up or Down button to increase/decrease the age. Press the SET button to select the displayed age.

8. Then step on the scale when it shows "0.0". The scale will automatically turn off after 60 seconds if no longer in use. Your personal data will be saved in the device's memory.


BODY WEIGHT MEASUREMENT

1. Place the scale on a flat, solid surface.
2. Turn on the scale by pressing the center of the platform.
3. Wait a few seconds and when "0.0" appears on the display, the scale is ready for weighing.
4. Step off the scale, the display will show the reading for a few more seconds. The scale automatically turns off.

MEASUREMENT OF FAT TISSUE CONTENT AND TOTAL BODY WATER.

Important information for users: The scale has 10 personal memory numbers (P1-P10). Before estimating body fat and water, you must first record your height, age, gender, normal/sport mode option and selecting the level of activity at a given number. The scale will then use factors to estimate the results. Once you have your personal information stored in memory, simply select. If the data changes, please enter it again. The accuracy of the results depends on where you stand on the scale. Take off your socks, get in position. Your feet must have maximum contact with the metal plate on the platform. Stay on the scale until the result appears on the display.


MEASUREMENT OF TISSUE FAT AND TOTAL BODY WATER USING USER PROFILE:

A. Selecting a user profile from memory:


1. Press the "SET" button to turn on the scale. The memory number (P1-P10) flashes.

2. Press the Up or Down buttons until your profile number appears. Wait until the screen shows "0.0".


3. When "0.0" appears on the screen, step on the scale with bare feet, placing your feet evenly on the scale platform with maximum contact with the metal electrodes.

4. Stand still while the scale measures your weight.

5. Your weight will be displayed for 2 seconds, then the scale will beep and start measuring your body composition. The screen will show a pattern of moving zeros "0000". Stand still on the scale.

6. After a few seconds, your BMI, body fat %, body water %, muscle mass %, KCAL and the estimated weight will be displayed.


7. The results will be repeated and the scale will turn off automatically.

B. Automatic user identification


After programming your profile in memory and weighing yourself, the scale will recognize you based on your weight. It does this based on the user's weight (it will detect all profiles within +/- 3kg/6.6lbs). If your weight differs from your previous one by +/- 3kg/6.6lbs, you will need to manually select a user profile.

1. Step onto the platform directly with your bare feet. Place your feet evenly, ensuring maximum contact with metal electrodes


2. Stand still while the scale measures your weight.

3. Your weight will be displayed for 2 seconds and then it will start measuring your body composition. During measurement, a moving zero pattern of "0000" will be visible on the screen. Stand still on the scale.

4. If there is only one user profile whose last scale reading is close to the person being weighed, the scale will automatically select that user's memory number. Your body composition measurement will then be displayed twice before the scale turns off.

5. If there is more than one user profile (e.g. 2 or 3) whose weight is similar to the weighed person (within 3 kg), the user profile numbers will appear on the screen. Step off the scale, press the "Up" or "Down" buttons to select the appropriate profile number.

HOW TO DELETE A USER PROFILE

1. Step on the scale. Your weight will appear.
2. The scale will try to select a user profile. If there are two or more conflicting profiles, the screen will display a group of conflicting profile numbers one by one within a few seconds.

3. Press the Up or Down buttons until the profile number you want to delete appears.

4. Press and hold the "SET" button for 3 seconds until the display shows "DEL"

5. The next time the scale is turned on, all information about the deleted user will be cleared.

BODY COMPOSITION:

Body mass index (BMI):

BMI is the ratio of body weight to height. The BMI table is shown below. A deviation above or below the BMI range may be considered a health risk.

BMI parameter (reference to WHO)				
	under weight	normal	over weight	obese
BMI	<18.5	18.5~24.9	25~29.9	>30

(ATTENTION! The table is an example. Each case is different and should be approached individually. If you do not have sufficient knowledge, consult a doctor or dietitian.)

Body fat %:

Body fat percentage represents your calculated total fat percentage

by adding essential fat and spare fat. Essential fats are necessary for various body functions.

Storage fat provides energy, protects against physical injury, and insulates the body to stay warm.

Amount of water in the body %:

The percentage of water in the body varies depending on gender: in healthy adult women it ranges from 45% to 60%, while in men it ranges from 50% to 65%. Slight fluctuations in body weight over a few hours are normal and result from water retention in the body.

Basal metabolic rate (Kcal):

Basal metabolic rate, BMR, is the minimum amount of calories the human body needs to survive, to perform basic functions such as circulation, breathing, digestion and maintaining body temperature balance.

SAFETY RULES:

1. Store and use the scale in a safe, dry place.
2. Do not disassemble the scale for purposes other than replacing the battery. .
3. After use, clean the scale with a soft, damp cloth.
4. Avoid excessive impact or vibration of the scale, such as dropping it on the floor.
5. This scale is designed to read weights from 5 to 225 kg (11 to 496 lb).
6. Don't store anything on the scale.
7. Use the scale in accordance with the recommendations contained in this manual.

QUESTIONS AND TROUBLESHOOTING

1. If the scale does not turn on or no weight is displayed?

Check if the scale has batteries and is turned on. If not, see the "Preparation" section before use"

2. If the scale displays "Lo"?

The batteries are exhausted. Replace all 3 used batteries.

3. If the scale is inaccurate or the weight changes?
Place your feet in the center of the platform and stand still. Wait until the weight stops before looking down on the screen. The floor may cause fluctuations. Use the scale on the hardest surface.

4. If the scale displays "88888" for more than 6 seconds?

Calibration error. Step off the scale and wait until the scale turns off automatically. Start weighing again by pressing firmly on the scale platform to recalibrate it. The display will show "0.0" and then turn off. The scale is ready for use again.

5. During measurement, the scale displays the message "Err"?

Overload warning. Remove the weight immediately; otherwise the device will be permanently damaged.

6. The scale displays "Err 1" when measuring body fat and water content?

Contact error. It is impossible to measure. Make sure you stand still on the scale, barefoot, and maintain maximum contact between your feet and the metal contacts. If not, see "MEASURING BODY FAT AND TOTAL BODY WATER." You may need to wet your feet to improve electrical contact.

7. The scale displays "Err L" when measuring body fat, water and muscle mass?

The estimated values are lower than the minimum values that the scale can estimate. Wetting your feet can help improve electrical contact.

8. The scale displays "Err H" when measuring body fat, water and muscle mass?

The estimated values are higher than the maximum values that the scale can estimate. Wetting your feet can help improve contact

9. What if I still can't solve the problem after taking the above steps?

If error messages still appear after following the steps above, remove and reinsert the batteries after 1 minute.

INFORMATION REGARDING BATTERY AND ACCUMULATION

- Keep batteries/rechargeable batteries away from children.
- Do not store batteries/rechargeable batteries in easily accessible places, as there is a risk that they may be swallowed by children or pets. If swallowed, consult a doctor as soon as possible.
- Leaky or damaged batteries/rechargeable batteries may cause burns when in contact with the skin, so always wear appropriate protective gloves when touching them.
- Batteries/rechargeable batteries must not be short-circuited, disassembled or thrown into fire. There is a risk of explosion!
- Ordinary disposable batteries should not be charged, as there is a risk of explosion! Only charge batteries intended for this purpose, using an appropriate charger.
- Charge the batteries before inserting them into the product
- Batteries should only be charged under the supervision of an adult.
- Do not mix alkaline, standard (carbon-zinc) and rechargeable (nickel-cadmium) batteries.
- Never mix new and old batteries.
- It is forbidden to interfere with the power supply cables - it may cause a short circuit.
- Exhausted batteries should be removed from the product.
- Do not throw used batteries into unsorted waste!
- Batteries should be disposed of in accordance with the regulations for the disposal of hazardous electronic waste.
- If batteries are not disposed of correctly, hazardous substances may pose a risk to human health or the environment. Recycling materials contributes to the conservation of natural resources.
- Always insert batteries/rechargeable batteries with polarity (plus/+ and minus/-).
- Type of batteries used in the product: 3XAAA

Watermarkly

Information for consumers on how to handle used electrical and electronic equipment:

Waste electrical and electronic equipment must not be disposed of with other household waste. Such equipment should be collected selectively, as indicated by the attached marking - crossed out wheeled waste container. Failure to comply with the above rule and improper disposal of used equipment may pose a threat to the natural environment and human health resulting from the presence of hazardous components in the equipment, such as: electrical wiring, plastic, switches, etc. To avoid the above threat, these components should be removed, collected and appropriately processed by specialized companies. Disassembling used equipment yourself is unacceptable. Households play an important role in contributing to reuse and recovery, including recycling of waste equipment. This is done in particular through active participation in the waste equipment collection system.

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