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Derma roller 3w1 0,25/0,5/1,5 mm

USER MANUAL

Derma roller 3in1 0.25 / 0.5 / 1.5 mm

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Derma-Walze 3in1 0,25 / 0,5 / 1,5 mm

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Derma váleček 3v1 0,25 / 0,5 / 1,5 mm

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Rouleau Derma 3 en 1 0,25 / 0,5 / 1,5 mm

INSTRUCCIONES DE OPERACION

Rodillo Derma 3 en 1 0,25 / 0,5 / 1,5 mm

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Rullo Derma 3in1 0,25 / 0,5 / 1,5 mm

MANUAL DE UTILIZARE

Rola Derma 3in1 0,25 / 0,5 / 1,5 mm

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Derma henger 3in1 0,25 / 0,5 / 1,5 mm

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Derma valček 3v1 0,25 / 0,5 / 1,5 mm

KASUTUSJUHEND

Derma rull 3in1 0,25 / 0,5 / 1,5 mm

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„Derma“ ritinėlis 3in1 0,25 / 0,5 / 1,5 mm

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Derma veltnis 3in1 0,25 / 0,5 / 1,5 mm

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EN

Ladies and Gentlemen, thank you for purchasing our product!
Before starting to use, please read the following instructions
for the correct use of the product. Please keep these instructions on
future and to follow its recommendations, as non-compliance may endanger life or health.

Important: KEEP FOR FUTURE!

- Please read it completely before use and follow all instructions in the manual.
- Do not use the product if any part is missing or damaged.
- Please keep this manual for future reference.

SPECIFICATION:

A home tool for performing micro-needle mesotherapy. Micro-needle mesotherapy is an innovative treatment that involves subjecting the skin to the deep action of revitalizing preparations. Their in-depth effect is obtained by introducing the preparations into the skin tissue with the help of micro-needles. Using the device evens out the skin tone, improves its firmness and smoothness, and reduces wrinkles. Roller comes with 3 different heads, from the largest to the smallest, for the body, face and under-eye area.

USE:

It is recommended to carry out a series of mesotherapy treatments at intervals of about 7-21 days. Before using the Derma Roller, the treatment site and the needle roller should be disinfected. Gently roll the selected part of the skin with a roller in different directions. Derma Roller can be gently pressed against the skin, but not too much discomfort. After the treatment, it is natural for slight, point

bleeding, reddening of the skin or swelling to occur. Before performing the next treatment, wait until the symptoms have completely disappeared.

STORAGE:

Store in sterile conditions. Keep away from children.

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