USER MANUAL

Sensorimotor pillow with a pump

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00022068 (MT-328-1+MT-328 SET)

EN

Ladies and Gentlemen, thank you for purchasing our product!
Before starting to use, please read the following instructions
for the correct use of the product. Please keep these instructions on
future and to follow its recommendations, as non-compliance may endanger life or health.

Important: KEEP FOR FUTURE REFERENCE!

- Please read it completely before use and follow all instructions in the manual.
- Do not use the product if any part is missing or damaged.
- Please keep this manual for future reference.

SPECIFICATION:

- Black colour.
- Diameter: 34 cm.

CHARACTERISTIC:

The pillow takes up little space, which makes it convenient to store and transport. It is an ideal solution for all people, regardless of age and gender. It can be used in fitness clubs, at home, in offices or in rehabilitation or physiotherapy centers. It has protrusions on both sides, so you can perform a variety of exercises depending on your needs and preferences.

INCLUDED:

- Sensory pillow.
- Pump with a needle.

CONTRAINDICATIONS:

- Pain during exercise.
- Conditions immediately after bone dislocation or fracture.
- Incomplete bone union.
- Inflammation of joints and periarticular tissues.
- Swelling.
- Fresh surgical scars or other wounds.

INDICATIONS:

- After removing the plaster cast.
- Labyrinth disorders in older people (frequent falls during everyday activities).
- Overload of stability, balance and motor coordination.
- Activation and strengthening of the postural muscles of the spine.

HOW TO USE:

- 1) Static sensorimotor training - the disc is used to perform coordination and balance exercises on an unstable surface and in open kinetic chains, e.g. maintaining balance by standing with one leg on a sensory pad, transferring body weight from leg to leg. It improves motor coordination in patients after injuries to the musculoskeletal system, increases the sense of stability in older people and reduces the risk of falls.
- 2) Activation of the postural muscles of the spine - sensory cushion during the so-called "dynamic sitting" forces the correct positioning of the back. Compared to standard sitting, the tension of postural muscles is normalized and the risk of spine overload is reduced.
- 3) Massage and relaxation of the feet - thanks to the sensory protrusions, the sensory receptors on the soles of the feet are irritated, their blood supply improves and, in response, the foot muscles relax.

CLEANING AND MAINTENANCE:

- 1) Properly inflate the pillow - there should be enough air in it to make it slightly convex. If you use it for exercise, inflate it so that when you stand on it you will have some trouble keeping your balance.
- 2) Clean the surface of the pillow regularly with a soft cloth or sponge, especially if used directly on the skin.
- 3) Avoid using chemicals such as solvents, detergents or cleaning agents as they may damage the surface of the pillow.
- 4) Avoid using aggressive chemicals such as solvents, detergents or cleaning agents as they may damage the surface of the pillow.
- 5) Store the pillow in a place with moderate temperature, avoid extreme heat and humidity conditions.

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