**soulima.**<sup>®</sup>

**USER MANUAL**

Ultrasonic peeling device

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5 902802 915917 >

EN

**Ladies and Gentlemen, thank you for purchasing our product!**  
**Before starting to use, please read the following instructions**  
**for the correct use of the product. Please keep these instructions on**  
**future and to follow its recommendations, as non-compliance may endanger life or health.**

**Important: KEEP FOR FUTURE REFERENCE!**

- Please read it completely before use and follow all instructions in the manual.
- Do not use the product if any part is missing or damaged.
- Please keep this manual for future reference.

**SPECIFICATION:**

- Black colour.
- Material: Plastic.
- USB cable; yes.
- Number of tips: 2.
- Dimensions: 17.5x5.5x1.5cm.
- Weight: 0.10kg.

**CHARACTERISTIC:**

The ultrasonic peeling device combines peeling and ultrasound - it cleanses the skin using saline, toner or water. Under the influence of ultrasonic vibrations, these bubbles burst, at the same time breaking up dead cells of the stratum corneum of the skin. The result is deep hydration, refreshment and alleviation of skin damage. Using the device, we can perform ultrasonic cavitation peeling at home.

**USER MANUAL:**

- 1) Remove the sealing film from the rejuvenating head before first use.
- 2) When the battery is low, the product will make a start-up sound, then the device needs to be charged.
- 3) To turn on the power, press and hold the power button for 2 seconds.
- 4) To activate the "rejuvenation" mode, press and hold the power button for 2 seconds. Then press once to select purification mode, press again to enter strong mode. Moisten the skin with water or lotion. Clean your face with the rejuvenating head
- 5) Press and hold the power button for 2 seconds to activate Moisturizing mode. Then press three times to choose humidification mode, press again to enter strong mode. Apply a cosmetic balm to the face (except the eyes) or apply a mask. Gently massage each part of your face for about 2-3 minutes. Then press five times to select Lifting mode, press again to enter Strong mode.

**SAFE USE RULES:**

- 1) Remove accessories such as necklace, earrings, bracelets, etc.
- 2) Please make sure that the power of the product is turned off before use.
- 3) Gently insert or pull out the charging cable to prevent damage to the port.
- 4) Store in a dry place.
- 5) To avoid damage to any parts, do not perform maintenance yourself.
- 6) When fully charged, unplug the charging cable.
- 7) Do not use or charge in high temperature or high humidity conditions.
- 8) Do not use the device with a damaged power cord.
- 9) Use the device with dry hands.
- 10) Do not press hard on the skin.
- 11) Do not use it in the same area for a long time.
- 12) Do not use around the eyes.
- 13) Do not use if you have any of the following symptoms: Facial neuropathy, Facial surgery or fracture with built-in surgery, cannot feel heat, allergies and sensitive skin, skin diseases and dermatitis, pregnant, congenital heart disease.

**CLEANING AND MAINTENANCE RULES:**

- 1) Cleaning should be systematic to prevent dust from accumulating and staying for a long time.
- 2) Clean with a dry, soft cloth. Dirt can be wiped with a damp cloth, and in the case of visible dirt (greasing) on the surface, a solution prepared from a few drops of ordinary dishwashing detergent dissolved in water can be used for cleaning, and then immediately dry the surfaces by wiping with a dry soft cloth. It is not recommended to soak with a lot of water.
- 3) Do not use cleaning agents containing abrasive substances, hard sponges, solvents or chemicals. Their use may cause damage and allergic reactions.
- 4) Do not store in a humid place such as a bathroom
- 5) Do not store in a place that is accessible to children

**Information for the consumer on the handling of waste electrical and electronic equipment:**

Waste electrical and electronic equipment may not be placed with other household waste. Such equipment should be collected selectively, which is indicated by the attached marking - a crossed-out wheeled waste container. Failure to comply with the above rule, in the case of improper disposal of used equipment, may pose a threat to the natural environment and human health, resulting from the presence of dangerous components in the equipment, such as: electrical wiring, plastic, switches, etc. To avoid the above threat, these components should be collected and properly processed by specialized companies. Self-disassembly of worn equipment is unacceptable. Households play an important role in contributing to reuse and recovery, including recycling of waste equipment. This is done in particular through active participation in the waste equipment collection system.

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